

CE EU Declaration of Conformity
 Notifying Medical Device Registration #202007742

CE

**Fiscal Year 2020
 CERTIFICATION OF REGISTRATION**

Model No.: ZST-B
 Size: For Adults 17.5X9.5cm
 Color: White Blue Pink

Shandong Zhushi Pharmaceutical Group Co., Ltd

Model No.: ZST-B
 Size: For Adults 17.5X9.5cm
 Color: White Blue Pink

Sun FDA

**Fiscal Year 2020
 CERTIFICATION OF REGISTRATION**

Model No.: ZST-B
 Size: For Adults 17.5X9.5cm
 Color: White Blue Pink

Model No.	Size	Color
20001	17.5x9.5cm	White
20002	17.5x9.5cm	Blue
20003	17.5x9.5cm	Pink

[Structure and composition]
 The product consists of surface layer, middle layer, bottom layer, mask belt and nose clip. The surface material is polypropylene spunbonded cloth, intermediate layer.

[Application]
 For clinical medical personnel to wear in the invasive operation process, covering the user's mouth, nose and jaw, in order to prevent pathogenic Direct transmission of microorganisms, body fluids, particulate matter, etc. Provides a physical screen.

[Attention]

1. Face mask can only be used once;
2. Replace the mask when it is wet;
3. Before entering the working area with the medical protective mask, the tightness should be checked;
4. If the mask is contaminated by the patient's blood and body fluids, it should be replaced in time;
5. Do not use after the package is damaged;
6. The product should be used as soon as possible after opening;
7. The product should be disposed according to the relevant provisions of medical waste after use.

SUNGO Europe B.V.
 VAT: NL857821659801
 Add: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands
 Contact: SUNGO Secretary
 Tel/Fax: +31 (0) 2021 11106
 E-mail: ec.repsungogroup.com

Shandong Zhushi Pharmaceutical Group Co., Ltd
 Add: No.6 Shande Road, Shan County, Heze City, Shandong, China 274300
 E-mail: zhushu-pharm.com

Made in China

FACE MASK

3ply with Earloop / Tie

CE ISO 10PCS

Shandong Zhushi Pharmaceutical Group Co., Ltd

INSTRUCTION

[Product Name] Face Mask
[Specification] 17.5x9.5cm
[Structure and composition]
 The product consists of surface layer, middle layer, bottom layer, mask belt and nose clip. The surface material is polypropylene spunbonded cloth, intermediate layer. The material is polypropylene spunbonded cloth, the bottom material is polypropylene spunbonded cloth, and the mask belt is polyester thread. The nose clip is made of polypropylene which can be bent and shaped.

[Application]
 For clinical medical personnel to wear in the invasive operation process, covering the user's mouth, nose and jaw, in order to prevent pathogenic Direct transmission of microorganisms, body fluids, particulate matter, etc. Provides a physical screen.

[Attention]

1. Face mask can only be used once;
2. Replace the mask when it is wet;
3. Before entering the working area with the medical protective mask, the tightness should be checked;
4. If the mask is contaminated by the patient's blood and body fluids, it should be replaced in time;
5. Do not use after the package is damaged;
6. The product should be used as soon as possible after opening;

[Storage conditions and shelf life]
 Shelf life of unopened product is 2 years from date of production when stored within temperature range of 20°C to 40°C and at less than 90% relative humidity.

[Manufacturer Name] Shandong Zhushi Pharmaceutical Group Co., LTD
[Address] No.6 Shande Road, Shan County, Heze City, Shandong, China.

CERTIFICATION

SHANDONG ZHUSHI PHARMACEUTICAL GROUP CO., LTD	
CERTIFICATION	
Product Name	Face Mask
Batch No.	LOT 20031902
Manufacture Date	19th Mar, 2020
Expiry Date	18th Mar, 2022
Inspector:	[Signature]
Period of Validity	2 Year