



LABSUN SELECTION

Mini Centrifuge

Mini Centrifuge

MicroCen 10

This mini centrifuge is ideal for quick spins, micro-filtration of samples, cell separations, and other routine laboratory procedures.




Features

- Ergonomically designed with a small footprint
- High grade plastic cover for clear visibility
- Smooth and quiet operation
- Automatic stop when lid is opened
- Digital control delivers precise speed
- Fits a variety of tubes with included accessories

Technical Parameters

Model	MicroCen 10
Product Name	Mini Centrifuge
Speed	10,000rpm
G-Force(RCF)	5,000xg
Noise Level	≤60dB
Body Color	Gray
Capacity	6×1.5ml / 6×0.5ml / 16×0.2ml
Accessories included in standard packaging	Rotor head: MLXF7(6×1.5/2.0ml) and MLXF3(16×0.2 ml) Adaptors: MLXF1(6×0.5ml) and MLXF2(6×0.2 ml)
Ambient Temperature	5°C~40°C
Operational Humidity	≤80%
Power Supply	AC110~240V, 50/60Hz
Max. Power	30W
Net Weight	≤1kg
Dimensions W × D × H	155×140×115mm
Safety Feature	Stops When Lid is Opened
Protection Level	IP21



Labsun GmbH

Am Sandfeld 13a
76149 Karlsruhe, Germany

Web: www.lab-sun.de

Email: info@lab-sun.de

Tel: +49 (0)721 47080525

Fax: +49 (0)721 47080537