



LABSUN SELECTION

Mini Vortex Mixer

Mini Vortex Mixer

VM-03UW

The mini vortex mixer is the ideal vortex mixer for every lab with its compact size, ease of use, and its economic price point.

The 3,000 rpm fixed speed ensures proper mixing and accurate results.




Features

- Powerful vortexing for tubes up to 30mm
- Extremely compact and easy operation
- Vortex head is made of durable material to ensure long life
- Refined steel base is designed to ensure product stability
- Choice of different colored working head

Technical Parameters

Model	VM-03UW
Product Name	Mini Vortex Mixer
Cap Color	Blue
Speed	3000 rpm
Speed Accuracy	±10%
Oscillation	Circular Oscillation
Orbital Diameter	∅ 4.5mm
Capacity	Max 80mL
Max Tube Diameter	30mm
Operating Environment	Ambient Temperature:5°C ~ 40°C Relative Humidity:≤80%
Power Supply	AC100-240V, 50/60Hz
Max.Power	5W
Dimensions Diameter x H	∅ 96 x 85 mm
Net Weight	0.6 kg
Protection Level	IP21



Labsun GmbH

Am Sandfeld 13a
76149 Karlsruhe, Germany

Web: www.lab-sun.de

Email: info@lab-sun.de

Tel: +49 (0)721 47080525

Fax: +49 (0)721 47080537